



## Intervertebral Disc Dehydration

- Degenerated intervertebral discs accumulate lactic acid causing pH levels to fall as low as 5.75
- Lactic acid alters cellular activity by down-regulating Proteoglycan (PG) synthesis and up-regulating enzymes that degrade the extracellular matrix.
- A low pH impairs the ability of PG to bind water causing the disc to degenerate and appear black under imaging



**Figure 1.** Image depicting the chemical formula of FehrFix hydrogel compound. Stable, dry material, made of cross-integrated hydro gel: Copolymers are based on methyl methacrylate and N-Vinylpyrrolidone



**Figure 2.** Image depicting FehrFix implants. Preoperatively implants are in a dry, non-inflated state and holds that rigid form and shape. Once implanted into the disc nucleus, they absorb natural body fluids as they hydrate and expand continuously to a controlled and predefined size and shape. Implants can absorb as much as 10 times their pre-op weight

## What is FehrFix?

- A non-surgical treatment for discogenic pain in the early stages of 'Degenerative Disc Disease' or when surgical intervention is undesirable or contraindicated
- FehrFix introduces a hydrophilic material (**hydrogel**) into the disc nucleus
- In-Situ volume expansion of the hydrogel restores natural hydration, pressure and stabilizes the pH in the nucleus to healthier levels
- FehrFix hydrates through the absorption of the body's own fluids
- The self-swelling hydrogel pins act upon the osmotic principle to gain volume

## FehrFix Hydrogel Implants

- High Biocompatibility – non-toxic
- Low Tissue Adherence
- Permeable to water soluble compounds
- Elasticity controlled by water
- Strength & Elasticity similar to tissues such as nucleus pulposus & cartilage

## Implant Procedure

- Non-surgical intervention under local anaesthesia
- Minimal trauma using cannulated needle
- Implant Inserted through a 17-gauge needle
- No pressure peaks occur
- ~15 minute procedure time

## FehrFix Clinical Information

### Study #1

- 13 Patients
- All with back and leg pain
- 10/13 (**77%**) Pain free or significant improvement
- 3 Patients required surgery (two after 6 months and one after 4 months)

### Study #2

- 17 Patients
- 16 with back and leg pain plus 1 with leg pain
- 15/17 (**88%**) Pain free or significant improvement
- 2/17 patients required a fusion due to disc collapse



### Indications

- Symptomatic lumbar DDD up to two segments.
- Dominant lumbar back pain with or without pseudo-radicular symptoms for more than 12 weeks, changes in MRI
- No improvement after minimum 12 week's conservative treatment with physiotherapy and pain drugs.
- No age limit
- Visible disc changes on MRI
- Minimal facet degeneration
- Annulus should be largely intact, HIZ (high intensity zone) is not contraindicated
- VAS Pain Score >3

### Contraindications

- Radiculopathy by nerve compression
- Stenosis, Spondylosis, Spondylolisthesis, acute fracture or Spondylitis at the treated segment.
- Modic type >2 or serious changes on endplates.
- Pregnancy